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No.MIICT0054EAC

1. CONFORMANCE STATEMENT OVERVIEW

Table 1-1 provides an overview of the network services.

Table 1-1
NETWORK SERVICES

| SOP Classes | User of Service (SCU) | Provider of Service (SCP) |
|---|-----------------------------|-----------------------------|
| Transfer | | |
| CT Image Storage | Yes | Yes* ¹ (COT-30D) |
| Secondary Capture Image Storage | Yes | Yes* ¹ (COT-30D) |
| Standalone Curve Storage | Yes | Yes* ¹ (COT-30D) |
| Enhanced CT Image Storage | Yes* ² (COT-45A) | Yes* ² (COT-45A) |
| Grayscale Softcopy Presentation State Storage | Yes | Yes* ¹ (COT-30D) |
| Enhanced SR Storage | Yes | Yes* ¹ (COT-30D) |
| X-Ray Radiation Dose SR Storage | Yes | Yes* ¹ (COT-30D) |
| Storage Commitment | | |
| Storage Commitment Push Model | Yes* ¹ (COT-41D) | No |
| Query/Retrieve | | |
| Patient Root Q/R Information Model – Find | No | Yes* ¹ (COT-34D) |
| Patient Root Q/R Information Model – Move | No | Yes* ¹ (COT-34D) |
| Study Root Q/R Information Model – Find | Yes* ¹ (COT-35D) | Yes* ¹ (COT-34D) |
| Study Root Q/R Information Model – Move | Yes* ¹ (COT-35D) | Yes* ¹ (COT-34D) |
| Patient/Study Only Information Model – Find | No | Yes* ¹ (COT-34D) |
| Patient/Study Only Information Model – Move | No | Yes* ¹ (COT-34D) |
| Workflow Management | | |
| Modality Worklist Information Model – Find | Yes* ¹ (COT-32D) | No |
| Modality Performed Procedure Step | Yes* ¹ (COT-33D) | No |
| Modality Performed Procedure Step Retrieve | Yes* ¹ (COT-33D) | No |

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